

50 Shades of Fluorescence to Follow Immune Response (Saffir) with spectral cytometry.

Introduction

Authors

Dimitri Popoff* Marc Barad* Juliette Desfrançois** **Cyrille Mionnet***

*Centre d'Immunologie de Marseille Luminy Marseille FRANCE CNRS UMR 7280, Inserm U1104, AMU-UM2. ** Cytek Biosciences

Methods

Whooping cough persists as an endemic disease (40million cases and 300000 child death per year worldwide) even in vaccinated populations. This disease is caused by the gram(-) bacterium Bordetella pertussis. This bacteria expresses multiple virulence factors acting in concert to facilitate its adherence, survival and proliferation in human respiratory tract. More, it has been shown, in a mouse model, that CD4⁺ Trm are generated during the time course of infection¹. Our objectives is to decipher the cellular mechanisms underlying the "cross-talk" between antigen presenting cells and T cells that allow the differentiation, maintenance and function of Trm at the tissue level. To do this aim we developped a 48 markers panel to follow the behavior of several immune populations in the mean time in the same tube.

• Mice were infected with a wild strain of Bordetella pertussis (Tahoma I). 15 days later they were sacrificed. After intra-cardiac perfusion with PBS, collected lungs were digested with a mix of Collagenase IV/DNAse I for 30 mn at 37°C. Lung cells were seeded at 3x10⁶ cells/well and stained in PBS, 2mM EDTA, 0,5% BSA, 40% brilliant violet stain buffer(BD Biosciences) and 10% Fc block at 37°C for 45mn. Cells washed 3 times before staining in PBS + live/dead marker (Zombie NIR) for 15mn. After 3 washes cells were resuspended in PBS for analysis on Spectral Cytometer Aurora (Cytek Biosciences).

Data analysis were performed using OMIQ.

ANTIGEN	CLONE	ANTIGEN	CLONE	ANTIGEN	CLONE	ANTIGEN	CLONE	ANTIGEN	CLONE
CD3	17A2	CD25	PC61	CD69	H1.2F3	CD196	29-2L17	Ly6G	IA8
CD4	GK1.5	CD26	H194-112	CD88	20/70	B220	RA3-6B2	MHC-II	M5/114.15.2
CD5	FABI 15U	CD27	LG3A10	CD90-2	30H12	BST2	927	NKI.I	PK136
CD8a	53-6.7	CD44	IM7	CD103	2 E7	CX3CRI	SAOLIFII	NKp46	29A1.4
CD8b	53-5.8	CD45	30F11	CD115	afs98	ЕрСАМ	G8.8	PD-I	29FIA12
CDIIb	MI/70	CD49a	HA31/8	CD138	281-2	F4/80	BM8	Siglec F	E50-2440
CDIIc	N418	CD49b	HMa2	CD154	MRI	IgM	RMM-I	XCRI	ZET
CD19	1D8	CD49d	PS/2	CD172a	P84	IgD	11-26c.2a	γδ TCR	GL3
CD23	B3B4	CD62L	MEL-14	CD183	CXCR3-183	Kirgl	2FI	Live/ Dead	Zombie NIR
CD24	M1/69	CD64	X54-5/7.1	CD192	SA203G11	Ly6C	HK1.4	AutoFluo	

Results CD3 CD4 PC24 PC25 PC26 PC27





Lower Panel:

Using supervised analysis (2 by 2) from CD45+->CD19⁻CMH-II⁻ -> CD90.2+CMH-II⁻ -> CD5+CD3+-> CD3+ $\gamma\delta T^-$ ->CD4+ or CD8a+ we finally are able to define a population of CD4+CD44+CD69+CD103+ (Trm) corresponding to 3,5% of total CD45⁺ cells. Idem CD8⁺ Trm represent 0,15% of CD45⁺ cells

Conclusion

Using Saffir we were able to follow the behavior of the majority of immune cells in response to Bordetella pertussis infection. Here, we highlighted the fact that using supervised or unsupervised analysis leaded to similar results in the frequency of both CD4+ and CD8+ Trm in lung of infected mice.

Reference: 1- Wilk, M. M. et al. Lung CD4 Tissue-Resident Memory T Cells Mediate Adaptive Immunity Induced by Previous Infection of Mice with Bordetella pertussis. J. Immunol. Baltim. Md 1950 199, 233–243 (2017)







